

Description and Use

Clinic is an antibacterial laundry powder designed for washing where temperatures required for thermal disinfection cannot be achieved, or are too high for the fabric type.

Clinic is low foaming so can be used in either front or top loading washing machines.

Features and Benefits

Clinic combines a high level of detergency with antibacterial action. When used as directed at 40°C, Clinic complies with the requirements of 'AS 4146:2000 Laundry practice' giving a plate count of <math><1 \text{ cfu/cm}^2</math>

Clinic contains activated oxygen bleach plus enzymes to give high level bleaching and stain removal properties. It also contains high levels of water softening agents to give excellent performance in all water conditions.

How to Use

For front loading washing machines use Clinic at a rate of 50g per 5kg capacity. This corresponds to 1 scoop or ¼ cup. Add extra if soiling is heavy.

In top loading machines use from 1.5 scoops in a 5kg machine to 2 scoops in a 7kg capacity machine.

Clinic has been tested at 40°C and gives excellent cleaning and antimicrobial performance at this temperature.

Do not add chlorine bleach to Clinic, or use for washing woollens or silk garments.

Technical Data

Clinic is an all in one laundry powder containing water softening agents, alkaline builders and biodegradable anionic and non-ionic detergents. It contains a protease enzyme to help remove protein based stains, plus an activated oxygen bleach system for stain removal and antimicrobial properties.

Clinic has a pH of 10.7 at a 1% solution.

Clinic has been independently tested at Silliker Microtech an independent TGA and NATA registered laboratory.

Tests were carried out in front and top loading washing machines using cloths impregnated with high levels of the following bacteria:

Staphylococcus aureus
Escherichia coli
Proteus vulgaris

The bacteria levels were measured before washing, after washing and after washing and drying. Even after just washing the bacteria were not detected.

Copies of the report are available on request

Pack Sizes

15kg bag	Code: 25466
4kg bucket	Code: 10396
1.5kg refill bag	Code: 87260

Dominant Colour Code



Hazard Information

Classified as hazardous according to Safe Work Australia criteria – R36 Irritating to eyes

Poisons Schedule:	Not classified
Dangerous Goods:	Not classified

Safety

Handle carefully to avoid creating dusts. Gloves are recommended when handling material. Keep out of reach of children

